# Rekommenderade metoder

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3. Beekmann SE, Diekema DJ, Chapin KC, Doern GV. Effects of rapid detection of bloodstream infections on length of hospitalization and hospital charges. *J Clin Microbiol.* 2003;41(7):3119-25.
4. Binkhamis K, Forward K. Effect of the initial specimen diversion technique on blood culture contamination rates. *J Clin Microbiol.* 2014;52(3):980-1.
5. Bouza E, Burillo A, Muñoz P, Guinea J, Marín M, Rodríguez-Créixems M. Mixed bloodstream infections involving bacteria and *Candida* spp. *The Journal of antimicrobial chemotherapy.* 2013;68(8):1881-8.
6. Cascone VJ, Cohen RS, Dodson NP, Cannon CM. Implications of culture collection after the first antimicrobial dose in septic emergency department patients. *Am J Emerg Med.* 2019;37(5):947-51.
7. Cattoir L, Claessens J, Cartuyvels R, Van den Abeele AM. How to achieve accurate blood culture volumes: the BD BACTEC FX blood volume monitoring system as a measuring instrument and educational tool. *European journal of clinical microbiology & infectious diseases : official publication of the European Society of Clinical Microbiology.* 2018;37(9):1621-6.
8. Chen IH, Nicolau DP, Kuti JL. Effect of Clinically Meaningful Antibiotic Concentrations on Recovery of *Escherichia coli* and *Klebsiella pneumoniae* Isolates from Anaerobic Blood Culture Bottles with and without Antibiotic Binding Resins. *J Clin Microbiol.* 2019;57(12).
9. Chen IH, Nicolau DP, Kuti JL. Recovery of Gram-Negative Bacteria from Aerobic Blood Culture Bottles Containing Antibiotic Binding Resins after Exposure to beta-Lactam and Fluoroquinolone Concentrations. *J Clin Microbiol.* 2019;57(10).
10. Cheng MP, Stenstrom R, Paquette K, Stabler SN, Akhter M, Davidson AC, et al. Blood Culture Results Before and After Antimicrobial Administration in Patients With Severe Manifestations of Sepsis: A Diagnostic Study. *Ann Intern Med.* 2019;171(8):547-54.
11. Chung Y, Kim IH, Han M, Kim HS, Kim HS, Song W, Kim JS. A comparative evaluation of BACT/ALERT FA PLUS and FN PLUS blood culture bottles and BD BACTEC Plus Aerobic and Anaerobic blood culture bottles for antimicrobial neutralization. *European journal of clinical microbiology & infectious diseases : official publication of the European Society of Clinical Microbiology.* 2019;38(12):2229-33.
12. Collin J. Svårigheter i att detektera *Candida albicans* i blododling vid samtidig bakteriemi. Examensarbete BMA062. Sahlgrenska akademien. 2019.
13. Dawson S. Blood culture contaminants. *J Hosp Infect.* 2014;87(1):1-10.
14. Dwivedi S, Bhalla R, Hoover DR, Weinstein MP. Discarding the initial aliquot of blood does not reduce contamination rates in intravenous-catheter-drawn blood cultures. *J Clin Microbiol.* 2009;47(9):2950-1.
15. Evans L, Rhodes A, Alhazzani W, Antonelli M, Coopersmith CM, French C, et al. Surviving Sepsis Campaign: International Guidelines for Management of Sepsis and Septic Shock 2021. *Critical care medicine.* 2021;49(11):e1063-e143.
16. Garcia RA, Spitzer ED, Beaudry J, Beck C, Diblasi R, Gilleeny-Blabac M, et al. Multidisciplinary team review of best practices for collection and handling of blood cultures to determine effective interventions for increasing the yield of true-positive bacteremias, reducing contamination, and eliminating false-positive central line-associated bloodstream infections. *Am J Infect Control.* 2015;43(11):1222-37.
17. Henning C, Aygül N, Dinnétz P, Wallgren K, Özenci V. Detailed Analysis of the Characteristics of Sample Volume in Blood Culture Bottles. *J Clin Microbiol.* 2019;57(8).
18. Klingspor L, Muhammed SA, Ozenci V. Comparison of the two blood culture systems, Bactec 9240 and BacT/Alert 3D, in the detection of *Candida* spp. and bacteria with

- polymicrobial sepsis. *European journal of clinical microbiology & infectious diseases* : official publication of the European Society of Clinical Microbiology. 2012;31(11):2983-7.
19. Lalezari A, Cohen MJ, Svinik O, Tel-Zur O, Sinvani S, Al-Dayem YA, et al. A simplified blood culture sampling protocol for reducing contamination and costs: a randomized controlled trial. *Clin Microbiol Infect*. 2020;26(4):470-4.
  20. Lamy B, Dargere S, Arendrup MC, Parienti JJ, Tattevin P. How to Optimize the Use of Blood Cultures for the Diagnosis of Bloodstream Infections? A State-of-the Art. *Front Microbiol*. 2016;7:697.
  21. Lamy B, Roy P, Carret G, Flandrois JP, Delignette-Muller ML. What is the relevance of obtaining multiple blood samples for culture? A comprehensive model to optimize the strategy for diagnosing bacteremia. *Clinical infectious diseases* : an official publication of the Infectious Diseases Society of America. 2002;35(7):842-50.
  22. Leyssene D, Gardes S, Vilquin P, Flandrois JP, Carret G, Lamy B. Species-driven interpretation guidelines in case of a single-sampling strategy for blood culture. *European journal of clinical microbiology & infectious diseases* : official publication of the European Society of Clinical Microbiology. 2011;30(12):1537-41.
  23. Meessen NEL, Jacobs JA. Blood volume in BACTECPLUS/F culture bottles sampled using the direct-draw technique. *Clinical Microbiology and Infection*. 1998;4(8):471-2.
  24. Menchinelli G, Liotti FM, Fiori B, De Angelis G, D'Inzeo T, Giordano L, et al. In vitro Evaluation of BACT/ALERT(R) VIRTUO(R), BACT/ALERT 3D(R), and BACTEC FX Automated Blood Culture Systems for Detection of Microbial Pathogens Using Simulated Human Blood Samples. *Front Microbiol*. 2019;10:221.
  25. Patton RG, Schmitt T. Innovation for reducing blood culture contamination: initial specimen diversion technique. *J Clin Microbiol*. 2010;48(12):4501-3.
  26. Péan de Ponfilly G, Benmansour H, Manda V, Lecorche E, Mougari F, Munier AL, et al. Impact of 24/7 loading of blood culture bottles in a new automated incubator on the diagnosis of bloodstream infections. *European journal of clinical microbiology & infectious diseases* : official publication of the European Society of Clinical Microbiology. 2021;40(12):2639-43.
  27. Pien BC, Mirrett S, Crews BR, Reller LB, Woods CW. Controlled clinical comparison of BacT/ALERT standard aerobic and standard anaerobic blood culture bottles inoculated directly or after transport in sodium polyanethol sulfonate tubes. *J Clin Microbiol*. 2007;45(4):1357-9.
  28. Rand KH, Beal SG, Rivera K, Allen B, Payton T, Lipori GP. Hourly Effect of Pretreatment With IV Antibiotics on Blood Culture Positivity Rate in Emergency Department Patients. *Open Forum Infect Dis*. 2019;6(5):ofz179.
  29. Riedel S. The value of postmortem microbiology cultures. *J Clin Microbiol*. 2014;52(4):1028-33.
  30. Ronnberg C, Mildh M, Ullberg M, Ozenci V. Transport time for blood culture bottles: underlying factors and its consequences. *Diagnostic microbiology and infectious disease*. 2013;76(3):286-90.
  31. Roth A, Wiklund AE, Pålsson AS, Melander EZ, Wullt M, Cronqvist J, et al. Reducing blood culture contamination by a simple informational intervention. *J Clin Microbiol*. 2010;48(12):4552-8.
  32. Rupp ME, Cavalieri RJ, Marolf C, Lyden E. Reduction in Blood Culture Contamination Through Use of Initial Specimen Diversion Device. *Clinical infectious diseases* : an official publication of the Infectious Diseases Society of America. 2017;65(2):201-5.

33. Sautter RL, Bills AR, Lang DL, Ruschell G, Heiter BJ, Bourbeau PP. Effects of delayed-entry conditions on the recovery and detection of microorganisms from BacT/ALERT and BACTEC blood culture bottles. *J Clin Microbiol.* 2006;44(4):1245-9.
34. Scheer CS, Fuchs C, Grundling M, Vollmer M, Bast J, Bohnert JA, et al. Impact of antibiotic administration on blood culture positivity at the beginning of sepsis: a prospective clinical cohort study. *Clin Microbiol Infect.* 2019;25(3):326-31.
35. Silver H SA. A practical and efficacious method for obtaining significant postmortem blood cultures. *American journal of clinical pathology.* 1969;52(4):433-7.
36. Syed S, Liss DT, Costas CO, Atkinson JM. Diversion Principle Reduces Skin Flora Contamination Rates in a Community Hospital. *Archives of pathology & laboratory medicine.* 2020;144(2):215-20.
37. Taniguchi T, Tsuha S, Shiiki S, Narita M. High positivity of blood cultures obtained within two hours after shaking chills. *International journal of infectious diseases : IJID : official publication of the International Society for Infectious Diseases.* 2018;76:23-8.
38. Yu D, Larsson A, Parke Å, Unge C, Henning C, Sundén-Cullberg J, et al. Single-Sampling Strategy vs. Multi-Sampling Strategy for Blood Cultures in Sepsis: A Prospective Non-inferiority Study. *Front Microbiol.* 2020;11:1639.
39. Zimmerman FS, Karamah H, Ben-Chetrit E, Zalut T, Assous M, Levin PD. Modification of blood test draw order to reduce blood culture contamination: a randomized clinical trial. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America.* 2019.
40. Åkerlund A, Petropoulos A, Malmros K, Tängdén T, Giske CG. Blood culture diagnostics: a Nordic multicentre survey comparison of practices in clinical microbiology laboratories. *Clin Microbiol Infect.* 2022;28(5):731.e1-.e7.

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1. Baron EJ, Miller JM, Weinstein MP, Richter SS, Gilligan PH, Thomson RB, Jr.,

- et al. A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2013 recommendations by the Infectious Diseases Society of America (IDSA) and the American Society for Microbiology (ASM)(a). *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2013;57(4):e22-e121.
2. Connell TG, Rele M, Cowley D, BATTERY JP, Curtis N. How reliable is a negative blood culture result? Volume of blood submitted for culture in routine practice in a children's hospital. *Pediatrics*. 2007;119(5):891-6.
  3. Dien Bard J, McElvania TeKippe E. Diagnosis of Bloodstream Infections in Children. *J Clin Microbiol*. 2016;54(6):1418-24.
  4. Dunne WM, Jr., Tillman J, Havens PL. Assessing the need for anaerobic medium for the recovery of clinically significant blood culture isolates in children. *The Pediatric infectious disease journal*. 1994;13(3):203-6.
  5. Howie SR. Blood sample volumes in child health research: review of safe limits. *Bull World Health Organ*. 2011;89(1):46-53.
  6. Huber S, Hetzer B, Crazzolara R, Orth-Höller D. The correct blood volume for paediatric blood cultures: a conundrum? *Clin Microbiol Infect*. 2020;26(2):168-73.
  7. Isaacman DJ, Karasic RB, Reynolds EA, Kost SI. Effect of number of blood cultures and volume of blood on detection of bacteremia in children. *The Journal of pediatrics*. 1996;128(2):190-5.
  8. Kellogg JA, Ferrentino FL, Goodstein MH, Liss J, Shapiro SL, Bankert DA. Frequency of low level bacteremia in infants from birth to two months of age. *The Pediatric infectious disease journal*. 1997;16(4):381-5.
  9. Kellogg JA, Manzella JP, Bankert DA. Frequency of low-level bacteremia in children from birth to fifteen years of age. *J Clin Microbiol*. 2000;38(6):2181-5.
  10. Lancaster DP, Friedman DF, Chiotos K, Sullivan KV. Blood Volume Required for Detection of Low Levels and Ultralow Levels of Organisms Responsible for Neonatal Bacteremia by Use of Bactec Peds Plus/F, Plus Aerobic/F Medium, and the BD Bactec FX System: an In Vitro Study. *J Clin Microbiol*. 2015;53(11):3609-13.
  11. Läkemedelsverket. Neonatal sepsis – ny behandlingsrekommendation. *Information från Läkemedelsverket* 2013;24:15-25.
  12. Melhus Å. Bakteriologisk diagnostik i neonatalperioden, *Information från Läkemedelsverket* 2013;24:44-8.
  13. Miller JM, Binnicker MJ, Campbell S, Carroll KC, Chapin KC, Gilligan PH, et al. A Guide to Utilization of the Microbiology Laboratory for Diagnosis of Infectious Diseases: 2018 Update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2018;67(6):813-6.
  14. Noel GJ, Laufer DA, Edelson PJ. Anaerobic bacteremia in a neonatal intensive care unit: an eighteen-year experience. *The Pediatric infectious disease journal*. 1988;7(12):858-62.
  15. Sarkar S, Bhagat I, DeCristofaro JD, Wiswell TE, Spitzer AR. A study of the role of multiple site blood cultures in the evaluation of neonatal sepsis. *Journal of perinatology : official journal of the California Perinatal Association*. 2006;26(1):18-22.
  16. Schelonka RL, Chai MK, Yoder BA, Hensley D, Brockett RM, Ascher DP. Volume of blood required to detect common neonatal pathogens. *The Journal of pediatrics*. 1996;129(2):275-8.
  17. Zaidi AK, Knaut AL, Mirrett S, Reller LB. Value of routine anaerobic blood cultures for pediatric patients. *The Journal of pediatrics*. 1995;127(2):263-8.

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#### Referenser Analys:

1. Akerlund A, Sundqvist M, Hanberger H, Ahren C, Serrander L, Giske CG. [Response times can be shortened in microbiological diagnosis of sepsis. Better opening hours in laboratories and active advice enables faster therapy]. *Lakartidningen*. 2015;112.
2. Almuhayawi M, Altun O, Abdulmajeed AD, Ullberg M, Ozenci V. The Performance of the Four Anaerobic Blood Culture Bottles BacT/ALERT-FN, -FN Plus, BACTEC-Plus and -Lytic in Detection of Anaerobic Bacteria and Identification by Direct MALDI-TOF MS. *PloS one*. 2015;10(11):e0142398.
3. Altun O, Almuhayawi M, Ullberg M, Ozenci V. Clinical evaluation of the FilmArray blood culture identification panel in identification of bacteria and yeasts from positive blood culture bottles. *J Clin Microbiol*. 2013;51(12):4130-6.
4. Altun O, Botero-Kleiven S, Carlsson S, Ullberg M, Ozenci V. Rapid identification of bacteria from positive blood culture bottles by MALDI-TOF MS following short-term incubation on solid media. *J Med Microbiol*. 2015;64(11):1346-52.
5. Arbetsmiljöverket. AFS 2018:4. Smittrisker. Arbetsmiljöverkets föreskrifter och allmänna råd om smittrisker.

6. Arbetsmiljöverket. AFS 2021:1. Smittrisker. Arbetsmiljöverkets föreskrifter om ändring i Arbetsmiljöverkets föreskrifter (AFS 2018:4) och allmänna råd om smittrisker.
7. Banerjee R, Ozenci V, Patel R. Individualized Approaches Are Needed for Optimized Blood Cultures. *Clinical infectious diseases* : an official publication of the Infectious Diseases Society of America. 2016;63(10):1332-9.
8. Barenfanger J, Graham DR, Kolluri L, Sangwan G, Lawhorn J, Drake CA, et al. Decreased mortality associated with prompt Gram staining of blood cultures. *American journal of clinical pathology*. 2008;130(6):870-6.
9. Baron EJ, Scott JD, Tompkins LS. Prolonged incubation and extensive subculturing do not increase recovery of clinically significant microorganisms from standard automated blood cultures. *Clinical infectious diseases* : an official publication of the Infectious Diseases Society of America. 2005;41(11):1677-80.
10. Beekmann SE, Diekema DJ, Chapin KC, Doern GV. Effects of rapid detection of bloodstream infections on length of hospitalization and hospital charges. *J Clin Microbiol*. 2003;41(7):3119-25.
11. Cornaglia G, Courcol R, Herrmann JL. *European Manual of Clinical Microbiology* 1st ed: European Society for Clinical Microbiology and Infections Diseases (ESCMID), Basel; 2012.
12. Edelmann A, Pietzcker T, Wellinghausen N. Comparison of direct disk diffusion and standard microtitre broth dilution susceptibility testing of blood culture isolates. *J Med Microbiol*. 2007;56(Pt 2):202-7.
13. Ekwall-Larson A, Fröding I, Mert B, Åkerlund A, Özenci V. Analytical Performance and Potential Clinical Utility of EUCAST Rapid Antimicrobial Susceptibility Testing in Blood Cultures after Four Hours of Incubation. *Microbiol Spectr*. 2023;11(2):e0500122.
14. Folkhälsomyndigheten. Packa provet rätt. 2020. Available from: <https://www.folkhalsomyndigheten.se/>
3. International Organization for Standardization. Susceptibility testing of infectious agents and evaluation of performance of antimicrobial susceptibility test devices—Part 1: Broth micro-dilution reference method for testing the in vitro activity of antimicrobial agents against rapidly growing aerobic bacteria involved in infectious diseases. 20776-1. 2019.
15. Johnsson ATA, Wong AYW, Özenci V. The impact of delayed analysis of positive blood cultures on the performance of short-term culture followed by MALDI-TOF MS. *J Microbiol Methods*. 2020;177:106027.
16. Jonasson E, Matuschek E, Kahlmeter G. The EUCAST rapid disc diffusion method for antimicrobial susceptibility testing directly from positive blood culture bottles. *Journal of Antimicrobial Chemotherapy*. 2020.
17. Lamy B, Sundqvist M, Idelevich EA. Bloodstream infections - Standard and progress in pathogen diagnostics. *Clin Microbiol Infect*. 2020;26(2):142-50.
4. Leber A. L. (Ed.) *Clinical Microbiology Procedures Handbook*. 4th ed: American Society for Microbiology; 2016.
18. Lin JN, Lai CH, Chen YH, Chang LL, Lu PL, Tsai SS, et al. Characteristics and outcomes of polymicrobial bloodstream infections in the emergency department: A matched case-control study. *Academic emergency medicine* : official journal of the Society for Academic Emergency Medicine. 2010;17(10):1072-9.
19. Matuschek E, Brown D, Kahlmeter G. Development of the EUCAST disk diffusion antimicrobial susceptibility testing method and its implementation in routine microbiology laboratories. *Clin Microbiol Infect* 2014 Apr;20(4):O255-66.

20. Pavlaki M, Poulakou G, Drimousis P, Adamis G, Apostolidou E, Gatselis NK, et al. Polymicrobial bloodstream infections: Epidemiology and impact on mortality. *Journal of global antimicrobial resistance*. 2013;1(4):207-12.
21. Petti CA, Bhally HS, Weinstein MP, Joho K, Wakefield T, Reller LB, Carroll KC. Utility of extended blood culture incubation for isolation of *Haemophilus*, *Actinobacillus*, *Cardiobacterium*, *Eikenella*, and *Kingella* organisms: a retrospective multicenter evaluation. *Journal of clinical microbiology*. 2006;44(1):257-9.
22. Public Health England. (2019). Investigation of blood cultures (for organisms other than *Mycobacterium* species). UK Standards for Microbiology Investigations. B 37 Issue 8.2. Available from: <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>. London: National Health Service (NHS). p. 1-51.
23. Savage TJ, Rao S, Joerger J, Ozonoff A, McAdam AJ, Sandora TJ. Predictive Value of Direct Disk Diffusion Testing from Positive Blood Cultures in a Children's Hospital and Its Utility in Antimicrobial Stewardship. *J Clin Microbiol*. 2021;59(6).
24. The European Committee on Antimicrobial Susceptibility Testing. EUCAST Disk Diffusion Method for Antimicrobial Susceptibility Testing. Version 12.0 (January 2024). Available at: [www.eucast.org](http://www.eucast.org).
25. World Health Organization. *Laboratory biosafety manual*. Geneva, WHO. 2020. Available from: <https://www.who.int/csr/resources/publications/biosafety/Biosafety7.pdf?ua=1>.
26. Åkerlund A, Jonasson E, Matuschek E, Serrander L, Sundqvist M, Kahlmeter G. EUCAST rapid antimicrobial susceptibility testing (RAST) in blood cultures: validation in 55 European laboratories. *The Journal of antimicrobial chemotherapy*. 2020;75(11):3230-8.
27. Åkerlund A, Petropoulos A, Malmros K, Tängdén T, Giske CG. Blood culture diagnostics: a Nordic multicentre survey comparison of practices in clinical microbiology laboratories. *Clin Microbiol Infect*. 2022;28(5):731.e1-.e7.

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#### Referenser Postanalys:

1. Akerlund A, Sundqvist M, Hanberger H, Ahren C, Serrander L, Giske CG. [Response times can be shortened in microbiological diagnosis of sepsis. Better opening hours in laboratories and active advice enables faster therapy]. *Lakartidningen*. 2015;112.
2. Banerjee R, Ozenci V, Patel R. Individualized Approaches Are Needed for Optimized Blood Cultures. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2016;63(10):1332-9.

3. Fang H, Ohlsson AK, Ullberg M, Ozenci V. Evaluation of species-specific PCR, Bruker MS, VITEK MS and the VITEK 2 system for the identification of clinical *Enterococcus* isolates. *European journal of clinical microbiology & infectious diseases* : official publication of the European Society of Clinical Microbiology. 2012;31(11):3073-7.
4. Haimi-Cohen Y, Shafinoori S, Tucci V, Rubin LG. Use of incubation time to detection in BACTEC 9240 to distinguish coagulase-negative staphylococcal contamination from infection in pediatric blood cultures. *The Pediatric infectious disease journal*. 2003;22(11):968-74.
5. Humphries RM, Kircher S, Ferrell A, Krause KM, Malherbe R, Hsiung A, Burnham CA. The Continued Value of Disk Diffusion for Assessing Antimicrobial Susceptibility in Clinical Laboratories: Report from the Clinical and Laboratory Standards Institute Methods Development and Standardization Working Group. *J Clin Microbiol*. 2018;56(8).
6. Munson EL, Diekema DJ, Beekmann SE, Chapin KC, Doern GV. Detection and treatment of bloodstream infection: laboratory reporting and antimicrobial management. *J Clin Microbiol*. 2003;41(1):495-7.
7. Ozenci V, Tegmark-Wisell K, Lundberg C, Wretling B. Rapid culture and identification: a practical method for early preliminary laboratory diagnosis of sepsis. *Clin Microbiol Infect*. 2008;14(2):177-80.
8. Public Health England. (2019). Investigation of blood cultures (for organisms other than *Mycobacterium* species). UK Standards for Microbiology Investigations. B 37 Issue 8.2. Available from: <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>. London: National Health Service (NHS). p. 1-51.
9. Ruus C, Skovbjerg S, Magnusson T, Snygg-Martin U, Studahl M, Andersson L-M. Tidig infektionskonsult gav effect vid *Staphylococcus aureus*-bakteriemi. *Läkartidningen*. 2018;3:88-91.
10. Smittskyddsförordning (2004:255). Stockholm: Socialdepartementet. Hämtad från: <http://rkrattsbaser.gov.se/sfst?bet=2004:255>.
11. Smittskyddslag (2004:168). Stockholm: Socialdepartementet. Hämtad från: <http://rkrattsbaser.gov.se/sfst?bet=2004:168>.
12. Åkerlund A, Petropoulos A, Malmros K, Tängdén T, Giske CG. Blood culture diagnostics: a Nordic multicentre survey comparison of practices in clinical microbiology laboratories. *Clin Microbiol Infect*. 2022;28(5):731.e1-.e7.

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1. Banerjee R, Ozenci V, Patel R. Individualized Approaches Are Needed for Optimized Blood Cultures. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2016;63(10):1332-9.
2. Barenfanger J, Graham DR, Kolluri L, Sangwan G, Lawhorn J, Drake CA, et al. Decreased mortality associated with prompt Gram staining of blood cultures. *American journal of clinical pathology*. 2008;130(6):870-6.
3. Cattoir L, Claessens J, Cartuyvels R, Van den Abeele AM. How to achieve accurate blood culture volumes: the BD BACTEC FX blood volume monitoring system as a measuring instrument and educational tool. *European journal of clinical microbiology & infectious diseases : official publication of the European Society of Clinical Microbiology*. 2018;37(9):1621-6.
4. Lamy B, Dargere S, Arendrup MC, Parienti JJ, Tattevin P. How to Optimize the Use of Blood Cultures for the Diagnosis of Bloodstream Infections? A State-of-the Art. *Front Microbiol*. 2016;7:697.
5. Lamy B, Ferroni A, Henning C, Cattoen C, Laudat P. How to: accreditation of blood cultures' proceedings. A clinical microbiology approach for adding value to patient care. *Clinical Microbiology and Infection*. 2018;24(9):956-63.
6. Lamy B, Sundqvist M, Idelevich EA. Bloodstream infections - Standard and progress in pathogen diagnostics. *Clin Microbiol Infect*. 2020;26(2):142-50.
7. Péan de Ponfilly G, Benmansour H, Manda V, Lecorche E, Mougari F, Munier AL, et al. Impact of 24/7 loading of blood culture bottles in a new automated incubator on the diagnosis of bloodstream infections. *European journal of clinical microbiology & infectious diseases : official publication of the European Society of Clinical Microbiology*. 2021;40(12):2639-43.
8. Public Health England. (2019). Investigation of blood cultures (for organisms other than Mycobacterium species). UK Standards for Microbiology Investigations. B 37 Issue 8.2. Available from: <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>. London: National Health Service (NHS). p. 1-51.

9. Sautter RL, Bills AR, Lang DL, Ruschell G, Heiter BJ, Bourbeau PP. Effects of delayed-entry conditions on the recovery and detection of microorganisms from BacT/ALERT and BACTEC blood culture bottles. *J Clin Microbiol.* 2006;44(4):1245-9.
10. Venturelli C, Righi E, Borsari L, Aggazzotti G, Busani S, Mussini C, et al. Impact of Pre-Analytical Time on the Recovery of Pathogens from Blood Cultures: Results from a Large Retrospective Survey. *PloS one.* 2017;12(1):e0169466.
11. Åkerlund A, Petropoulos A, Malmros K, Tängdén T, Giske CG. Blood culture diagnostics: a Nordic multicentre survey comparison of practices in clinical microbiology laboratories. *Clin Microbiol Infect.* 2022;28(5):731.e1-.e7.

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1. Acuña M, O’Ryan M, Cofré J, Alvarez I, Benadof D, Rodríguez P, et al. Differential time to positivity and quantitative cultures for noninvasive diagnosis of catheter-related blood stream infection in children. *The Pediatric infectious disease journal*. 2008;27(8):681-5.
2. Afshari A, Schrenzel J, Ieven M, Harbarth S. Bench-to-bedside review: Rapid molecular diagnostics for bloodstream infection--a new frontier? *Crit Care*. 2012;16(3):222.
3. Armstrong CW, Mayhall CG, Miller KB, Newsome HH, Jr., Sugerman HJ, Dalton HP, et al. Clinical predictors of infection of central venous catheters used for total parenteral nutrition. *Infect Control Hosp Epidemiol*. 1990;11(2):71-8.
4. Ben-Ami R, Weinberger M, Orni-Wasserlauff R, Schwartz D, Itzhaki A, Lazarovitch T, et al. Time to blood culture positivity as a marker for catheter-related candidemia. *J Clin Microbiol*. 2008;46(7):2222-6.
5. Bong JJ, Kite P, Ammori BJ, Wilcox MH, McMahon MJ. The use of a rapid in situ test in the detection of central venous catheter-related bloodstream infection: a prospective study. *JPEN J Parenter Enteral Nutr*. 2003;27(2):146-50.
6. Bouza E, Alcalá L, Muñoz P, Martín-Rabadán P, Guembe M, Rodríguez-Creixems M. Can microbiologists help to assess catheter involvement in candidaemic patients before removal? *Clinical microbiology and infection : the official publication of the European Society of Clinical Microbiology and Infectious Diseases*. 2013;19(2):E129-35.
7. Bouza E, Alvarado N, Alcalá L, Sánchez-Conde M, Pérez MJ, Muñoz P, et al. A prospective, randomized, and comparative study of 3 different methods for the diagnosis of intravascular catheter colonization. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2005;40(8):1096-100.
8. Bouza E, Burillo A, Muñoz P. Catheter-related infections: diagnosis and intravascular treatment. *Clinical microbiology and infection : the official publication of the European Society of Clinical Microbiology and Infectious Diseases*. 2002;8(5):265-74.
9. Bouzidi H, Emirian A, Marty A, Chachaty E, Laplanche A, Gachot B, Blot F. Differential time to positivity of central and peripheral blood cultures is inaccurate for the diagnosis of *Staphylococcus aureus* long-term catheter-related sepsis. *J Hosp Infect*. 2018;99(2):192-9.
10. Brun-Buisson C, Abrouk F, Legrand P, Huet Y, Larabi S, Rapin M. Diagnosis of central venous catheter-related sepsis. Critical level of quantitative tip cultures. *Arch Intern Med*. 1987;147(5):873-7.
11. Catton JA, Dobbins BM, Kite P, Wood JM, Eastwood K, Sugden S, et al. In situ diagnosis of intravascular catheter-related bloodstream infection: a comparison of quantitative culture, differential time to positivity, and endoluminal brushing. *Critical care medicine*. 2005;33(4):787-91.

12. Chaiyakunapruk N, Veenstra DL, Lipsky BA, Sullivan SD, Saint S. Vascular catheter site care: the clinical and economic benefits of chlorhexidine gluconate compared with povidone iodine. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2003;37(6):764-71.
13. Erb S, Frei R, Schreggenberger K, Dangel M, Nogarth D, Widmer AF. Sonication for diagnosis of catheter-related infection is not better than traditional roll-plate culture: a prospective cohort study with 975 central venous catheters. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2014;59(4):541-4.
14. Falagas ME, Ierodiakonou V, Alexiou VG. Clinical practice of obtaining blood cultures from patients with a central venous catheter in place: an international survey. *Clinical microbiology and infection : the official publication of the European Society of Clinical Microbiology and Infectious Diseases*. 2009;15(7):683-6.
15. Fenollar F, Raoult D. Molecular diagnosis of bloodstream infections caused by non-cultivable bacteria. *Int J Antimicrob Agents*. 2007;30 Suppl 1:S7-15.
16. Fernández-Cruz A, Martín-Rabadán P, Suárez-Salas M, Rojas-Wettig L, Pérez MJ, Guinea J, et al. Is it feasible to diagnose catheter-related candidemia without catheter withdrawal? *Med Mycol*. 2014;52(5):491-7.
17. Gahlot R NC, Kumar V, Yadav G, Anupurba S. Catheter-related bloodstream infections. *Int J Crit Illn Inj Sci*. 2014;4(2):162-7.
28. Garrison RN, Wilson MA. Intravenous and central catheter infections. *Surg Clin North Am*. 1994;74(3):557-70.
19. Gominet M, Compain F, Beloin C, Lebeaux D. Central venous catheters and biofilms: where do we stand in 2017? *APMIS : acta pathologica, microbiologica, et immunologica Scandinavica*. 2017;125(4):365-75.
20. Gowardman JR, Jeffries P, Lassig-Smith M, Stuart J, Jarrett P, Deans R, et al. A comparative assessment of two conservative methods for the diagnosis of catheter-related infection in critically ill patients. *Intensive Care Med*. 2013;39(1):109-16.
21. Gristina AG, Costerton JW. Bacterial adherence to biomaterials and tissue. The significance of its role in clinical sepsis. *J Bone Joint Surg Am*. 1985;67(2):264-73.
22. Gueembe M, Martin-Rabadan P, Cruces R, Perez-Granda MJ, Bouza E. Corrigendum to "Sonicating multi-lumen sliced catheter tips after the roll-plate technique improves the detection of catheter colonization in adults" [*J. Microbiol. Methods* 122 (2016) 20-22]. *Journal of microbiological methods*. 2016;127:242.
23. Gueembe M, Martin-Rabadan P, Cruces R, Perez Granda MJ, Bouza E. Sonicating multi-lumen sliced catheter tips after the roll-plate technique improves the detection of catheter colonization in adults. *Journal of microbiological methods*. 2016;122:20-2.
24. Hamedi A, Mamoury GA, Akhlaghi F. A comparative study of blood culture sampling from umbilical catheter line versus peripheral site. *Acta medica Iranica*. 2010;48(4):231-3.
25. Hammarskjöld F, Berg S, Hanberger H, Taxbro K, Malmvall BE. Sustained low incidence of central venous catheter-related infections over six years in a Swedish hospital with an active central venous catheter team. *American journal of infection control*. 2014;42(2):122-8.
26. Hansen S SF, Behnke M, Carsaw H, Heczko P, Klavs I, Lyytikäinen O, Palomar M, Riesenfeld Orn I, Savey A, Szilagyi E, Valinteliene R, Fabry J, Gastmeier P. National influences on catheter-associated bloodstream infection rates: practices among national surveillance networks participating in the European HELICS project. *J Hosp Infect*. 2009;1(1):66-73.

27. Hoiby N, Bjarnsholt T, Moser C, Bassi GL, Coenye T, Donelli G, et al. ESCMID guideline for the diagnosis and treatment of biofilm infections 2014. *Clinical microbiology and infection : the official publication of the European Society of Clinical Microbiology and Infectious Diseases*. 2015;21 Suppl 1:S1-25.
28. Hollenbeak CS. The cost of catheter-related bloodstream infections: implications for the value of prevention. *J Infus Nurs*. 2011;34(5):309-13.
29. Kaasch AJ, Rieg S, Hellmich M, Kern WV, Seifert H. Differential time to positivity is not predictive for central line-related *Staphylococcus aureus* bloodstream infection in routine clinical care. *J Infect*. 2014;68(1):58-61.
30. Kumar A, Sharma RM, Jaideep CN, Hazra N. Diagnosis of central venous catheter-related bloodstream infection without catheter removal: A prospective observational study. *Med J Armed Forces India*. 2014;70(1):17-21.
31. Maki DG, Weise CE, Sarafin HW. A semiquantitative culture method for identifying intravenous-catheter-related infection. *N Engl J Med*. 1977;296(23):1305-9.
32. Monsen T, Lövgren E, Widerström M, Wallinder L. In Vitro Effect of Ultrasound on Bacteria and Suggested Protocol for Sonication and Diagnosis of Prosthetic Infections. *Journal of Clinical Microbiology*. 2009;47(8):2496-501.
33. Murray PR, Masur H. Current approaches to the diagnosis of bacterial and fungal bloodstream infections in the intensive care unit. *Critical care medicine*. 2012;40(12):3277-82.
34. National Nosocomial Infections Surveillance (NNIS) System Report, Data Summary from January 1990-May 1999, issued June 1999. A report from the NNIS System. *American journal of infection control*. 1999;27(6):520-32.
35. O'Grady NP, Alexander M, Burns LA, Dellinger EP, Garland J, Heard SO, et al. Guidelines for the prevention of intravascular catheter-related infections. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2011;52(9):e162-93.
36. Orihuela-Martin J, Rodriguez-Nunez O, Morata L, Cardozo C, Puerta-Alcalde P, Hernandez-Meneses M, et al. Performance of differential time to positivity as a routine diagnostic test for catheter-related bloodstream infections: a single-centre experience. *Clin Microbiol Infect*. 2020;26(3):383.e1-e7.
37. Parienti JJ, du Cheyron D, Ramakers M, Malbruny B, Leclercq R, Le Coutour X, Charbonneau P. Alcoholic povidone-iodine to prevent central venous catheter colonization: A randomized unit-crossover study. *Critical care medicine*. 2004;32(3):708-13.
38. Public Health England. (2017). Investigation of intravascular cannulae and associated specimens. UK Standards for Microbiology Investigations. B 20 Issue 6.1. Available from: <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>. London: National Health Service (NHS). p. 1-24
39. Quittnat Pelletier F, Joarder M, Poutanen SM, Lok CE. Evaluating Approaches for the Diagnosis of Hemodialysis Catheter-Related Bloodstream Infections. *Clin J Am Soc Nephrol*. 2016;11(5):847-54.
40. Raad I, Hanna H, Maki D. Intravascular catheter-related infections: advances in diagnosis, prevention, and management. *The Lancet Infectious diseases*. 2007;7(10):645-57.
41. Raad I, Hanna HA, Alakech B, Chatzinikolaou I, Johnson MM, Tarrand J. Differential time to positivity: a useful method for diagnosing catheter-related bloodstream infections. *Ann Intern Med*. 2004;140(1):18-25.
42. Sherertz RJ, Raad, II, Belani A, Koo LC, Rand KH, Pickett DL, et al. Three-year experience with sonicated vascular catheter cultures in a clinical microbiology laboratory. *Journal of clinical microbiology*. 1990;28(1):76-82.

43. Siegman-Igra Y, Anglim AM, Shapiro DE, Adal KA, Strain BA, Farr BM. Diagnosis of vascular catheter-related bloodstream infection: a meta-analysis. *Journal of clinical microbiology*. 1997;35(4):928-36.
44. Stempel JM, Farmakiotis D, Tarrand JJ, Kontoyiannis DP. Time-to-reporting of blood culture positivity and central venous catheter-associated *Candida glabrata* fungemia in cancer patients. *Diagnostic microbiology and infectious disease*. 2016;85(3):391-3.
45. Tighe MJ, Kite P, Thomas D, Fawley WN, McMahon MJ. Rapid diagnosis of catheter-related sepsis using the acridine orange leukocyte cytospin test and an endoluminal brush. *JPEN J Parenter Enteral Nutr*. 1996;20(3):215-8.
46. Towns ML, Jarvis WR, Hsueh PR. Guidelines on blood cultures. *J Microbiol Immunol Infect*. 2010;43(4):347-9.
47. von Baum H, Philippi P, Geiss HK. Acridine-orange leucocyte cytospin (AOLC) test as an in-situ method for the diagnosis of central venous catheter (CVC)-related sepsis in adult risk patients. *Zentralbl Bakteriol*. 1998;287(1-2):117-23.
48. Wu H, Moser C, Wang HZ, Hoiby N, Song ZJ. Strategies for combating bacterial biofilm infections. *Int J Oral Sci*. 2015;7(1):1-7.
49. Yousif A, Jamal MA, Raad I. Biofilm-based central line-associated bloodstream infections. *Adv Exp Med Biol*. 2015;830:157-79.

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1. Alestig K, Hogevis H, Olaison L. Infective endocarditis: a diagnostic and therapeutic challenge for the new millennium. *Scand J Infect Dis.* 2000;32(4):343-56.
2. Baddour LM, Epstein AE, Erickson CC, Knight BP, Levison ME, Lockhart PB, et al. Update on cardiovascular implantable electronic device infections and their management: a scientific statement from the American Heart Association. *Circulation.* 2010;121(3):458-77.
3. Baron EJ, Scott JD, Tompkins LS. Prolonged incubation and extensive subculturing do not increase recovery of clinically significant microorganisms from standard automated blood cultures. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America.* 2005;41(11):1677-80.
4. Cahill TJ, Prendergast BD. Infective endocarditis. *Lancet.* 2016;387(10021):882-93.
5. Calabrese F, Carturan E, Thiene G. Cardiac infections: focus on molecular diagnosis. *Cardiovasc Pathol.* 2010;19(3):171-82.
6. Cockerill FR, 3rd, Wilson JW, Vetter EA, Goodman KM, Torgerson CA, Harmsen WS, et al. Optimal testing parameters for blood cultures. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America.* 2004;38(12):1724-30.

7. Dayer MJ, Jones S, Prendergast B, Baddour LM, Lockhart PB, Thornhill MH. Incidence of infective endocarditis in England, 2000-13: a secular trend, interrupted time-series analysis. *Lancet*. 2015;385(9974):1219-28.
8. Durack DT, Lukes AS, Bright DK. New criteria for diagnosis of infective endocarditis: utilization of specific echocardiographic findings. Duke Endocarditis Service. *Am J Med*. 1994;96(3):200-9.
9. Fowler VG, Durack DT, Selton-Suty C, Athan E, Bayer AS, Chamis AL, et al. The 2023 Duke-International Society for Cardiovascular Infectious Diseases Criteria for Infective Endocarditis: Updating the Modified Duke Criteria. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2023;77(4):518-26.
10. Gould FK, Denning DW, Elliott TS, Foweraker J, Perry JD, Prendergast BD, et al. Guidelines for the diagnosis and antibiotic treatment of endocarditis in adults: a report of the Working Party of the British Society for Antimicrobial Chemotherapy. *The Journal of antimicrobial chemotherapy*. 2012;67(2):269-89.
11. Greub G, Lepidi H, Rovey C, Casalta JP, Habib G, Collard F, et al. Diagnosis of infectious endocarditis in patients undergoing valve surgery. *Am J Med*. 2005;118(3):230-8.
12. Hogevis H, Olaison L, Andersson R, Lindberg J, Alestig K. Epidemiologic aspects of infective endocarditis in an urban population. A 5-year prospective study. *Medicine (Baltimore)*. 1995;74(6):324-39.
13. Lee A, Mirrett S, Reller LB, Weinstein MP. Detection of bloodstream infections in adults: how many blood cultures are needed? *Journal of clinical microbiology*. 2007;45(11):3546-8.
14. Li J, Plorde JJ, Carlson LG. Effects of volume and periodicity on blood cultures. *Journal of clinical microbiology*. 1994;32(11):2829-31.
15. Li JS, Sexton DJ, Mick N, Nettles R, Fowler VG, Jr., Ryan T, et al. Proposed modifications to the Duke criteria for the diagnosis of infective endocarditis. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2000;30(4):633-8.
16. Liesman RM, Pritt BS, Maleszewski JJ, Patel R. Laboratory Diagnosis of Infective Endocarditis. *Journal of clinical microbiology*. 2017;55(9):2599-608.
17. Lobmaier IV, Vege A, Gaustad P, Rognum TO. Bacteriological investigation--significance of time lapse after death. *European journal of clinical microbiology & infectious diseases : official publication of the European Society of Clinical Microbiology*. 2009;28(10):1191-8.
18. Mermel LA, Maki DG. Detection of bacteremia in adults: consequences of culturing an inadequate volume of blood. *Ann Intern Med*. 1993;119(4):270-2.
19. Morris AJ, Drinkovic D, Pottumarthy S, Strickett MG, MacCulloch D, Lambie N, Kerr AR. Gram stain, culture, and histopathological examination findings for heart valves removed because of infective endocarditis. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2003;36(6):697-704.
20. Munoz P, Bouza E, Marin M, Alcalá L, Rodríguez Creixems M, Valerio M, Pinto A. Heart valves should not be routinely cultured. *Journal of clinical microbiology*. 2008;46(9):2897-901.
21. Ortiz C, Lopez J, Garcia H, Sevilla T, Revilla A, Vilacosta I, et al. Clinical classification and prognosis of isolated right-sided infective endocarditis. *Medicine (Baltimore)*. 2014;93(27):e137.
22. Petti CA, Bhally HS, Weinstein MP, Joho K, Wakefield T, Reller LB, Carroll KC. Utility of extended blood culture incubation for isolation of *Haemophilus*,

- Actinobacillus, Cardiobacterium, Eikenella, and Kingella organisms: a retrospective multicenter evaluation. *Journal of clinical microbiology*. 2006;44(1):257-9.
23. Public Health England. (2019). Investigation of blood cultures (for organisms other than Mycobacterium species). UK Standards for Microbiology Investigations. B 37 Issue 8.2. Available from: <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>. London: National Health Service (NHS). p. 1-51.
  24. Razimi R. MP. ABC om Infektiös endokardit. *Läkartidningen*. 2019;116:1-5.
  25. Röver C, Greub G, Lepidi H, Casalta JP, Habib G, Collart F, Raoult D. PCR detection of bacteria on cardiac valves of patients with treated bacterial endocarditis. *Journal of clinical microbiology*. 2005;43(1):163-7.
  26. Silver H SA. A practical and efficacious method for obtaining significant postmortem blood cultures. *American journal of clinical pathology*. 1969;52(4):433-7.
  27. Svenska infektionsläkarföreningen. Programgruppen för endokardit 2019–2021. Vårdprogram, Infektiös endokardit version 7, 2021. Available from: [http://infektion.net/wp-content/uploads/2017/05/vardprogram\\_infektios\\_endokardit.pdf](http://infektion.net/wp-content/uploads/2017/05/vardprogram_infektios_endokardit.pdf)
  28. Ternhag A, Cederstrom A, Torner A, Westling K. A nationwide cohort study of mortality risk and long-term prognosis in infective endocarditis in Sweden. *PloS one*. 2013;8(7):e67519.
  29. Vondracek M, Sartipy U, Aufwerber E, Julander I, Lindblom D, Westling K. 16S rDNA sequencing of valve tissue improves microbiological diagnosis in surgically treated patients with infective endocarditis. *The Journal of infection*. 2011;62(6):472-8.
  30. Werner M. Blood culture negative endocarditis. Göteborg: Göteborg University; 2006.
  31. Werner M, Andersson R, Olaison L, Hogevik H. A 10-year survey of blood culture negative endocarditis in Sweden: aminoglycoside therapy is important for survival. *Scand J Infect Dis*. 2008;40(4):279-85.
  32. Working Party of the British Society for Antimicrobial Chemotherapy. Guidelines for the diagnosis and antibiotic treatment of endocarditis in adults: a report of the Working Party of the British Society for Antimicrobial Chemotherapy. *J Antimicrob Chemother* 67: British Society for Antimicrobial Chemotherapy. 2012. p. 269-89.

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1. Banerjee R, Ozenci V, Patel R. Individualized Approaches Are Needed for Optimized Blood Cultures. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2016;63(10):1332-9.
2. Bell M, Bogar C, Plante J, Rasmussen K, Winters S. Effectiveness of a Novel Specimen Collection System in Reducing Blood Culture Contamination Rates. *J Emerg Nurs*. 2018;44(6):570-5.
3. Binkhamis K, Forward K. Effect of the initial specimen diversion technique on blood culture contamination rates. *J Clin Microbiol*. 2014;52(3):980-1.
4. Brown E, Wenzel RP, Hendley JO. Exploration of the microbial anatomy of normal human skin by using plasmid profiles of coagulase-negative staphylococci: search for the reservoir of resident skin flora. *J Infect Dis*. 1989;160(4):644-50.
5. Cascone VJ, Cohen RS, Dodson NP, Cannon CM. Implications of culture collection after the first antimicrobial dose in septic emergency department patients. *Am J Emerg Med*. 2019;37(5):947-51.
6. Cattoir L, Claessens J, Cartuyvels R, Van den Abeele AM. How to achieve accurate blood culture volumes: the BD BACTEC FX blood volume monitoring system as a measuring instrument and educational tool. *European journal of clinical microbiology & infectious diseases : official publication of the European Society of Clinical Microbiology*. 2018;37(9):1621-6.
7. Chen IH, Nicolau DP, Kuti JL. Effect of Clinically Meaningful Antibiotic Concentrations on Recovery of Escherichia coli and Klebsiella pneumoniae Isolates from Anaerobic Blood Culture Bottles with and without Antibiotic Binding Resins. *J Clin Microbiol*. 2019;57(12).
8. Chen IH, Nicolau DP, Kuti JL. Recovery of Gram-Negative Bacteria from Aerobic Blood Culture Bottles Containing Antibiotic Binding Resins after Exposure to beta-Lactam and Fluoroquinolone Concentrations. *J Clin Microbiol*. 2019;57(10).
9. Cheng MP, Stenstrom R, Paquette K, Stabler SN, Akhter M, Davidson AC, et al. Blood Culture Results Before and After Antimicrobial Administration in Patients With Severe Manifestations of Sepsis: A Diagnostic Study. *Ann Intern Med*. 2019;171(8):547-54.
10. Chung Y, Kim IH, Han M, Kim HS, Kim HS, Song W, Kim JS. A comparative evaluation of BACT/ALERT FA PLUS and FN PLUS blood culture bottles and BD

- BACTEC Plus Aerobic and Anaerobic blood culture bottles for antimicrobial neutralization. *European journal of clinical microbiology & infectious diseases* : official publication of the European Society of Clinical Microbiology. 2019;38(12):2229-33.
11. Dwivedi S, Bhalla R, Hoover DR, Weinstein MP. Discarding the initial aliquot of blood does not reduce contamination rates in intravenous-catheter-drawn blood cultures. *J Clin Microbiol.* 2009;47(9):2950-1.
  12. Garcia RA, Spitzer ED, Beaudry J, Beck C, Diblasi R, Gilleeny-Blabac M, et al. Multidisciplinary team review of best practices for collection and handling of blood cultures to determine effective interventions for increasing the yield of true-positive bacteremias, reducing contamination, and eliminating false-positive central line-associated bloodstream infections. *Am J Infect Control.* 2015;43(11):1222-37.
  13. Gonsalves WI, Cornish N, Moore M, Chen A, Varman M. Effects of volume and site of blood draw on blood culture results. *J Clin Microbiol.* 2009;47(11):3482-5.
  14. Grice EA, Kong HH, Renaud G, Young AC, Bouffard GG, Blakesley RW, et al. A diversity profile of the human skin microbiota. *Genome Res.* 2008;18(7):1043-50.
  15. Hall KK, Lyman JA. Updated review of blood culture contamination. *Clinical microbiology reviews.* 2006;19(4):788-802.
  16. Henning C, Aygül N, Dinnétz P, Wallgren K, Özenci V. Detailed Analysis of the Characteristics of Sample Volume in Blood Culture Bottles. *J Clin Microbiol.* 2019;57(8).
  17. Jonsson B, Nyberg A, Henning C. Theoretical aspects of detection of bacteraemia as a function of the volume of blood cultured. *Apmis.* 1993;101(8):595-601.
  18. Lalezari A, Cohen MJ, Svinik O, Tel-Zur O, Sinvani S, Al-Dayem YA, et al. A simplified blood culture sampling protocol for reducing contamination and costs: a randomized controlled trial. *Clin Microbiol Infect.* 2020;26(4):470-4.
  19. Lamy B, Dargere S, Arendrup MC, Parienti JJ, Tattevin P. How to Optimize the Use of Blood Cultures for the Diagnosis of Bloodstream Infections? A State-of-the Art. *Front Microbiol.* 2016;7:697.
  20. Lamy B, Roy P, Carret G, Flandrois JP, Delignette-Muller ML. What is the relevance of obtaining multiple blood samples for culture? A comprehensive model to optimize the strategy for diagnosing bacteremia. *Clinical infectious diseases* : an official publication of the Infectious Diseases Society of America. 2002;35(7):842-50.
  21. Leysse D, Gardes S, Vilquin P, Flandrois JP, Carret G, Lamy B. Species-driven interpretation guidelines in case of a single-sampling strategy for blood culture. *European journal of clinical microbiology & infectious diseases* : official publication of the European Society of Clinical Microbiology. 2011;30(12):1537-41.
  22. Li J, Plorde JJ, Carlson LG. Effects of volume and periodicity on blood cultures. *Journal of clinical microbiology.* 1994;32(11):2829-31.
  23. Liumbruno GM, Catalano L, Piccinini V, Pupella S, Grazzini G. Reduction of the risk of bacterial contamination of blood components through diversion of the first part of the donation of blood and blood components. *Blood Transfus.* 2009;7(2):86-93.
  24. Lockhart PB, Brennan MT, Sasser HC, Fox PC, Paster BJ, Bahrani-Mougeot FK. Bacteremia associated with toothbrushing and dental extraction. *Circulation.* 2008;117(24):3118-25.
  25. Meessen NEL, Jacobs JA. Blood volume in BACTECPLUS/F culture bottles sampled using the direct-draw technique. *Clinical Microbiology and Infection.* 1998;4(8):471-2.
  26. Menchinelli G, Liotti FM, Fiori B, De Angelis G, D'Inzeo T, Giordano L, et al. In vitro Evaluation of BACT/ALERT(R) VIRTUO(R), BACT/ALERT 3D(R), and BACTEC FX Automated Blood Culture Systems for Detection of Microbial Pathogens Using Simulated Human Blood Samples. *Front Microbiol.* 2019;10:221.

27. Mitteregger D, Barousch W, Nehr M, Kundi M, Zeitlinger M, Makristathis A, Hirschl AM. Neutralization of antimicrobial substances in new BacT/Alert FA and FN Plus blood culture bottles. *J Clin Microbiol.* 2013;51(5):1534-40.
28. Palarasah Y, Skjoedt MO, Vitved L, Andersen TE, Skjoedt K, Koch C. Sodium polyanethole sulfonate as an inhibitor of activation of complement function in blood culture systems. *J Clin Microbiol.* 2010;48(3):908-14.
29. Patton RG, Schmitt T. Innovation for reducing blood culture contamination: initial specimen diversion technique. *J Clin Microbiol.* 2010;48(12):4501-3.
30. Pien BC, Mirrett S, Crews BR, Reller LB, Woods CW. Controlled clinical comparison of BacT/ALERT standard aerobic and standard anaerobic blood culture bottles inoculated directly or after transport in sodium polyanethol sulfonate tubes. *J Clin Microbiol.* 2007;45(4):1357-9.
31. Rand KH, Beal SG, Rivera K, Allen B, Payton T, Lipori GP. Hourly Effect of Pretreatment With IV Antibiotics on Blood Culture Positivity Rate in Emergency Department Patients. *Open Forum Infect Dis.* 2019;6(5):ofz179.
32. Roth A, Wiklund AE, Pålsson AS, Melander EZ, Wullt M, Cronqvist J, et al. Reducing blood culture contamination by a simple informational intervention. *J Clin Microbiol.* 2010;48(12):4552-8.
33. Rupp ME, Cavalieri RJ, Marolf C, Lyden E. Reduction in Blood Culture Contamination Through Use of Initial Specimen Diversion Device. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America.* 2017;65(2):201-5.
34. Syed S, Liss DT, Costas CO, Atkinson JM. Diversion Principle Reduces Skin Flora Contamination Rates in a Community Hospital. *Archives of pathology & laboratory medicine.* 2020;144(2):215-20.
35. Vincent JL, Sakr Y, Sprung CL, Ranieri VM, Reinhart K, Gerlach H, et al. Sepsis in European intensive care units: results of the SOAP study. *Critical care medicine.* 2006;34(2):344-53.
36. Yu D, Larsson A, Parke Å, Unge C, Henning C, Sundén-Cullberg J, et al. Single-Sampling Strategy vs. Multi-Sampling Strategy for Blood Cultures in Sepsis: A Prospective Non-inferiority Study. *Front Microbiol.* 2020;11:1639.
37. Zimmerman FS, Karamah H, Ben-Chetrit E, Zalut T, Assous M, Levin PD. Modification of blood test draw order to reduce blood culture contamination: a randomized clinical trial. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America.* 2019.

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patient care. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. 2013;56(9):1284-92.

5. Diekema D, Arbefeville S, Boyken L, Kroeger J, Pfaller M. The changing epidemiology of healthcare-associated candidemia over three decades. *Diagnostic microbiology and infectious disease*. 2012;73(1):45-8.
6. Ericson EL, Klingspor L, Ullberg M, Ozenci V. Clinical comparison of the Bactec Mycosis IC/F, BacT/Alert FA, and BacT/Alert FN blood culture vials for the detection of candidemia. *Diagnostic microbiology and infectious disease*. 2012;73(2):153-6.
7. Hockey LJ, Fujita NK, Gibson TR, Rotrosen D, Montgomerie JZ, Edwards JE, Jr. Detection of fungemia obscured by concomitant bacteremia: in vitro and in vivo studies. *J Clin Microbiol*. 1982;16(6):1080-5.
8. Klingspor L. Svampinfektioner (invasiva), diagnostik. Stockholm: Internetmedicin.se, [uppdaterad 2018-01-21; citerad 2019-09-04]. Available from: <https://www.internetmedicin.se/page.aspx?id=615>
9. Klingspor L, Lindback E, Ullberg M, Ozenci V. Seven years of clinical experience with the Yeast Traffic Light PNA FISH: Assay performance and possible implications on antifungal therapy. *Mycoses*. 2018;61(3):179-85.
10. Klingspor L, Muhammed SA, Ozenci V. Comparison of the two blood culture systems, Bactec 9240 and BacT/Alert 3D, in the detection of *Candida* spp. and bacteria with polymicrobial sepsis. *European journal of clinical microbiology & infectious diseases: official publication of the European Society of Clinical Microbiology*. 2012;31(11):2983-7.
11. Klingspor L, Ullberg M, Rydberg J, Kondori N, Serrander L, Swanberg J, et al. Epidemiology of fungaemia in Sweden: A nationwide retrospective observational survey. *Mycoses*. 2018;61(10):777-85.
12. Klotz SA, Chasin BS, Powell B, Gaur NK, Lipke PN. Polymicrobial bloodstream infections involving *Candida* species: analysis of patients and review of the literature. *Diagnostic microbiology and infectious disease*. 2007;59(4):401-6.
13. Mellinghoff SC, Cornely OA, Jung N. Essentials in *Candida* bloodstream infection. *Infection*. 2018;46(6):897-9.
14. Pappas PG, Lionakis MS, Arendrup MC, Ostrosky-Zeichner L, Kullberg BJ. Invasive candidiasis. *Nature reviews Disease primers*. 2018;4:18026.
15. Public Health England. (2019). Investigation of blood cultures (for organisms other than *Mycobacterium* species). UK Standards for Microbiology Investigations. B 37 Issue 8.2. Available from: <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>. London: National Health Service (NHS). p. 1-51.
16. Referensgruppen för Antimykotika. Laboratoriediagnostik och farmakologisk monitorering av invasiva svampinfektioner i Sverige 2022. Available from: <https://infektion.net/wp-content/uploads/2021/08/ram-laboratoriediagnostik-och-farmakologisk-monitorering-av-invasiva-svampinfektioner-i-sverige.pdf>.
17. Sjölin J HL. Svampinfektioner (invasiva), behandling [Internet]. Göteborg Internetmedicin.se, [uppdaterad 2018-01-14; citerad 2019-08-30]. Available from: <https://www.internetmedicin.se/page.aspx?id=614>.
18. Stempel JM, Farmakiotis D, Tarrand JJ, Kontoyiannis DP. Time-to-reporting of blood culture positivity and central venous catheter-associated *Candida glabrata* fungemia in cancer patients. *Diagnostic microbiology and infectious disease*. 2016;85(3):391-3.

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